

SILVERSG™

Antibacterial Wound Dressing Gel

The clear solution in wound care™

SILVERSG™ Antibacterial Wound Dressing Gel contains silver that in laboratory tests has been shown to inhibit the growth of microorganisms such as:

- *Staphylococcus aureus*
- *Pseudomonas aeruginosa*
- *Escherichia coli*
- Antibiotic-resistant bacteria
- MRSA and VRE
- Fungi such as *Candida albicans*.

SILVERSG™ is a water-based gel indicated for use in the management of:

- First and second degree burns
- Wounds such as stasis ulcers, pressure ulcers, diabetic ulcers
- Lacerations
- Abrasions
- Skin tears
- Surgical incision sites
- Device insertion site wounds
- Graft sites
- Donor sites

SILVERSG™ is a clear gel that does not obscure the wound from view and there is no discoloration.

0.5 OZ (14.174g) - REF 10290

1 OZ (28.349g) - REF 10291

3 OZ (85.05g) - REF 10292

Caution: Federal law restricts this product to sale by or on the order of a licensed healthcare practitioner.



SteriGear®

www.sterigear.com

800-398-3259 - Provo, UT